

Methods

The system holds over 150 user definable applications. IL offers a comprehensive menu that includes reagents, calibrators and controls for a wide range of validated applications covering routine chemistry, specific proteins, rheumatology... Available dedicated kits are listed in the tables below.

ROUTINE CHEMISTRY

Test	Method	Cat. No.	Kit Configuration	Tests/Kit	Calibrators	Formulation
Acid Phosphatase (Total & Non Prostatic)	α -naphthylphosphate	00018482300	5x6 (NP) + 5x6 (T) ml	200	ReferrIL G	Liquid
Albumin	Bromocresol Green	00018481800	8x10 mL	250	ReferrIL G	Liquid
Alkaline Phosphatase	IFCC	00018480600	16x6.5 mL	400	Factor	Liquid
ALT/GPT	IFCC	00018480700	6x26 + 6x6.5 mL	696	Factor	Dry
ALT/GPT	IFCC	00018257440	8x48 + 4x24 mL	1714	Factor	Liquid
Amylase	CNP-Maltotrioidase	00018481500	12x8 mL	320	Factor	Liquid
AST/GOT	IFCC	00018480800	6x26 + 6x6.5 mL	696	Factor	Dry
AST/GOT	IFCC	00018257540	8x48 + 4x24 mL	1714	Factor	Liquid
Calcium	o-cresolphthalein	00018482800	6x30 + 6x30 mL	600	ReferrIL G	Liquid
Cholesterol	Chol. Esterase Chol. Oxidase / Trinder	00018480200	6x35 mL	656	ReferrIL G	Liquid
Cholinesterase	BTC	00018482000	6x10 + 4x30 mL	363	Factor	Dry
CK	IFCC	00018482100	8x10 + 4x5 mL	320	Factor	Liquid
CK-MB	Immunoinhibition	00018482200	6x10 + 6x2.5 mL	250	Factor	Liquid
Creatinine	Jaffè	00018480900	6x16 + 6x16 mL	640	ReferrIL G	Liquid
Direct Bilirubin	Modified Jendrassik-Grof	00018481100	4x33 + 8x10 mL	471	ReferrIL G	Liquid
γ -GT	IFCC	00018480300	4x17 mL	226	Factor	Dry
γ -GT	Szasz / IFCC	00018257640	6x40 + 6x10 mL	1000	Factor	Liquid
Glucose	Glucose Oxidase/Trinder	00018480000	6x35 mL	700	ReferrIL G	Liquid
HDL Cholesterol	Direct Homogeneous	00018481400	4x30 + 4x10 mL	400	ReferrIL HDL	Liquid
LDL Cholesterol	Direct Homogeneous	00018482700	6x9 + 6x3 mL	180	ReferrIL LDL	Liquid
Iron	Ferene	00018481700	4x30 + 4x2.5 mL	400	ReferrIL G	Liquid
Magnesium	Xilidyl Blue	00018481600	8x10 mL	228	ReferrIL G	Liquid
Phosphorus	Phosphomolybdate Complex	00018481900	6x10 + 6x10 mL	400	ReferrIL G	Liquid
Total Bilirubin	Modified Jendrassik-Grof	00018481000	4x33 + 8x10 mL	440	ReferrIL G	Liquid
Total Protein	Biuret	00018481300	6x30 mL	600	ReferrIL G	Liquid
Triglycerides	GPO/Trinder	00018480500	6x30 mL	600	ReferrIL G	Liquid
Urea	Urease/GLDH	00018480400	8x26 + 8x6.5 mL	866	ReferrIL G	Liquid
Uric Acid	Uricase/Trinder	00018480100	6x33 mL	660	ReferrIL G	Dry

IMMUNOTURBIDIMETRY

Test	Method	Cat. No.	Kit Configuration	Tests/Kit	Calibrators	Formulation
Quantex CRP	Latex-immunoturbidimetry	T3000-2321	4x13 + 4x6.5 mL	162	CRP Standard	Liquid
Quantex ASO	Latex-immunoturbidimetry	T3000-2319	4x13 + 4x2 mL	160	ASO Standard	Liquid
Quantex RF	Latex-immunoturbidimetry	T3000-2320	4x13 + 4x2 mL	160	RF Standard	Liquid
Quantex Ferritin	Latex-immunoturbidimetry	T3000-2328	8x4 + 2x4 mL	129	Ferritin standard	Liquid
Quantex IgA	Immunoturbidimetric	T3000-2322	1x30 + 4x6 mL	96	Protein standard	Liquid
Quantex IgG	Immunoturbidimetric	T3000-2323	1x30 + 4x6 mL	96	Protein standard	Liquid
Quantex IgM	Immunoturbidimetric	T3000-2324	1x30 + 4x5 mL	100	Protein standard	Liquid
Quantex Transferrin	Immunoturbidimetric	T3000-2325	2x21 + 2x2 mL	133	Protein standard	Liquid

CALIBRATORS

	Description	Cat. No.	Kit Configuration	Parameters	Composition
Routine Chemistry	ReferrIL G	00018257000	10x3 mL	Substrates and Electrolytes (no ISE)	Lyoph. bovine based
	ReferrIL HDL Cholesterol	00018255800	2x3 mL	HDL Chol	Lyoph. human based
	ReferrIL LDL Cholesterol	00018256100	2x1 mL	LDL Chol	Lyoph. human based
	ISE CALIBRATOR A	00018469700	10x100 mL	Na, K, Cl	Acqueous Standard
	ISE CALIBRATOR B	00018469800	4x20 mL	Na, K, Cl	Acqueous Standard
Immunochemistry	Quantex RF Plus standard	T3000-2252	6x1 mL	RF	Liquid human based
	Quantex Ferritin standard multipoint	T3000-2223	5x1 mL	Ferritin	Liquid human based
	Quantex CRP Plus standard	T3000-2093	6x1 mL	CRP	Liquid human based
	Quantex ASO Plus standard	T3000-2312	6x1 mL	ASO	Liquid human based
	Quantex Protein standard multipoint	T3000-2128	5x1 mL	IgA, IgG, IgM, alpha1-AGP, Transferrin	Liquid human based

CONTROLS

	Description	Cat. No.	Kit Configuration	Parameters	Composition
Routine Chemistry	SeraChem Control Level 1	00018162412	12x5 mL	Routine Chemistry	Lyoph. human based
	SeraChem Control Level 2	00018162512	12x5 mL	Routine Chemistry	Lyoph. human based
Immunochemistry	Quantex ASO-CRP-RF control I	T3000-2069	6x1 mL	ASO, CRP, RF	Liquid human based
	Quantex ASO-CRP-RF control II	T3000-2070	3x1 mL	ASO, CRP, RF	Liquid human based
	Quantex Ferritin/Myoglobin/IgE control I-II	T3000-2222	2x3 mL	Ferritin	Liquid human based
	Quantex Protein control I-II	T3000-2122	3x1 + 3x1 mL	IgA, IgG, IgM, alpha1-AGP, Transferrin	Liquid human based

The ILab Test Menu also includes: Glucose Hexokinase, Lipase, LD-LDH, Soluble Transferrin, Immunoglobulin E, Myoglobin, Complement 3, Complement 4, Microalbumin, Haptoglobin, HbA1c. For detailed information, please, refer to ILab 300 Plus Application Manual P/N 00018461800.