

ROUTINE CHEMISTRY

Test	Method	Cat. No.	Kit Configuration	Tests/Kit	Calibrators	Formulation
Acid Phosphatase (Total & Non Prostatic)	α -naphthylphosphate	00018482300	5x6(7)+5x8(NP) mL	284	Factor	Dry
Albumin	Bromocresol Green	00018250040	5x95 mL	1938	ReferriL G	Liquid
Albumin	Bromocresol Green	00018481800	8x10 mL	326	ReferriL G	Liquid
Alkaline Phosphatase	IFCC	00018252140	4x60 mL	1411	Factor	Liquid
Alkaline Phosphatase	IFCC	00018480600	16x7.5 mL	705	Factor	Liquid
ALT/GPT	IFCC	00018257440	8x49+4x24 mL	2742	Factor	Liquid
ALT/GPT	IFCC	00018252240	10x20+10x11 mL	2750	Factor	Dry
Amylase	CNP-Maltotrioxide	00018255340	10x20 mL	1000	Factor	Liquid
Amylase	CNP-Maltotrioxide	00018481500	12x8 mL	480	Factor	Liquid
AST/GOT	IFCC	00018257540	8x49+4x24 mL	2742	Factor	Liquid
AST/GOT	IFCC	00018252440	10x20+10x11 mL	2750	Factor	Dry
Calcium	o-cresolphthalein	00018250340	5x100+5x100 mL	2631	ReferriL G	Liquid
Cholesterol	Chol. Esterase Chol. Oxidase / Trinder	00018250540	5x50 mL	2500	ReferriL G	Liquid
Cholinesterase	BTC	00018254340	4x90+2x90 mL	2769	Factor	Dry
CK	IFCC	00018257340	6x24+3x12 mL	800	Factor	Liquid
CK	IFCC	00018482100	8x10+4x5 mL	444	Factor	Liquid
CK-MB	Immunoinhibition	00018482200	6x10+6x2.5 mL	333	Factor	Liquid
Creatinine	Jaffe	00018257240	10x13+10x24 mL	2626	ReferriL G	Liquid
Direct Bilirubin	Modified Jendrassik-Grof	00018254740	4x90+4x75.75 mL	3030	ReferriL G	Liquid
γ -GT	IFCC	00018254240	4x90 mL	2482	Factor	Dry
γ -GT	Szasz / IFCC	00018257640	6x40 + 6x10 mL	1714	Factor	Liquid
Glucose Hexokinase	Glucose Hexokinase	00018250740	8x90 mL	3428	ReferriL G	Liquid
Glucose Oxidase	Glucose Oxidase/Trinder	00018250840	5x95 mL	3166	ReferriL G	Liquid
HDL Cholesterol	Direct Homogeneous	00018255740	2x60+2x20 mL	444	ReferriL HDL	Liquid
Iron	Ferene	00018255040	4x90+2x15 mL	1440	ReferriL G	Liquid
LDH-P	Pyruvate-lactate	00018254940	4x90+4x20 mL	1777	Factor	Liquid
LDL Cholesterol	Direct Homogeneous	00018256040	4x18+4x6 mL	266	ReferriL LDL	Liquid
Lipase	Colorimetric, kinetic	00018254440	6x10+2x10+1x3 mL	307	Lipase calibrator (kit included)	Liquid
Magnesium	Xilidyl Blue	00018256240	5x20 mL	285	ReferriL G	Liquid
Phosphorus	Phosphomolybdate Complex	00018251240	5x45+5x45 mL	1500	ReferriL G	Liquid
Phosphorus	Phosphomolybdate Complex	00018481900	6x10+6x10 mL	400	ReferriL G	Liquid
Total Bilirubin	Modified Jendrassik-Grof	00018254640	4x100+4x75.75 mL	2885	ReferriL G	Liquid
Total Protein	Buret	00018251440	5x100 mL	2360	ReferriL G	Liquid
Triglycerides	GPO/Trinder	00018255640	5x100 mL	1851	ReferriL G	Liquid
Urea	Urease/GLDH	00018255440	8x100 mL	2962	ReferriL G	Liquid
Uric Acid	Uricase/Trinder	00018251740	5x90 mL	2500	ReferriL G	Dry

IMMUNODIAGNOSTIC

Test	Method	Cat. No.	Kit Configuration	Tests/Kit	Calibrators	Formulation
Quantex α 1 Acid Glycoprotein	Immunoturbidimetric	T3000-2179	1x100+1x4 mL	200	Protein standard	Liquid
Quantex α 1 Antitrypsin	Immunoturbidimetric	T3000-2306	1x35+2x4 mL	160	Protein standard	Liquid
Quantex ASO	Latex-immunoturbidimetry	T3000-2304	1x90+2x6 mL	240	ASO Standard	Liquid
Quantex β 2 Microglobulin	Immunoturbidimetric	T3000-2307	1x40+4x3 mL	100	β -2 Microglobulin standard	Liquid
Quantex C3	Immunoturbidimetric	T3000-2308	1x60+4x12 mL	252	Protein standard	Liquid
Quantex C4	Immunoturbidimetric	T3000-2309	1x60+4x12 mL	252	Protein standard	Liquid
Quantex CRP	Latex-immunoturbidimetry	T3000-2209	1x80+4x10 mL	250	CRP Standard	Liquid
Quantex CRP ultra sensitive	Latex-immunoturbidimetry	T3000-2297	1x13+1x13 mL	123	CRPus Standard	Liquid
Quantex D-DIMER	Immunoturbidimetric	T3000-8501	2x13+2x3 mL	100	D-DIMER standard	Dry
Quantex Ferritin	Latex-immunoturbidimetry	T3000-2271	1x100+3x7 mL	381	Ferritin standard	Liquid
Quantex Haptoglobin	Immunoturbidimetric	T3000-2305	1x35+2x3.5 mL	155	Protein standard	Liquid
Quantex HbA1c	Immunoturbidimetric	T3000-2314	1x50+1x28+1x14+1x14 mL	143	HbA1c standard	Liquid
Quantex IgA	Immunoturbidimetric	T3000-2268	1x90+4x16 mL	241	Protein standard	Liquid
Quantex IgE	Latex-immunoturbidimetry	T3000-2238	2x16+2x5 mL	117	IgE standard	Liquid
Quantex IgG	Immunoturbidimetric	T3000-2269	1x90+4x16 mL	241	Protein standard	Liquid
Quantex IgM	Immunoturbidimetric	T3000-2270	1x90+4x16 mL	304	Protein standard	Liquid
Quantex Lp(a)	Latex-immunoturbidimetry	T3000-2292	2x11+2x4mL	123	Lp(a) standard	Liquid
Quantex Microalbumin	Immunoturbidimetric	T3000-2272	1x75+1x4 mL	200	Microalbumin standard	Liquid
Quantex Myoglobin	Latex-immunoturbidimetry	T3000-2329	2x12+2x6 mL	126	Myoglobin standard	Liquid
Quantex RF	Latex-immunoturbidimetry	T3000-2255	1x75+2x6 mL	300	RF Standard	Liquid
Quantex Rubella	Latex-immunoturbidimetry	T3000-2180	2x5+2x4 mL	76	Rubella standard	Dry
Quantex sTfR	Latex-immunoturbidimetry	T3000-2300	1x45+1x4 mL	133	sTfR standard	Dry
Quantex Toxo	Latex-immunoturbidimetry	T3000-2181	2x14+2x3.5 mL	116	Toxo standard	Liquid
Quantex Transferrin	Immunoturbidimetric	T3000-2172	2x100+2x4 mL	400	Protein standard	Liquid

THERAPEUTIC DRUG MONITORING

Test	Method	Cat. No.	Kit Configuration	Tests/Kit	Calibrators	Formulation
Quantex Amikacin	Latex-immunoturbidimetry	T3000-2317	2x12.5+2x3.8 mL	94	Amikacin standard	Liquid
Quantex Carbamazepine	Latex-immunoturbidimetry	T3000-2277	2x11+2x3 mL	93	Carbamazepine standard	Liquid
Quantex Digitoxin	Latex-immunoturbidimetry	T3000-2276	2x10.5+2x6.5 mL	182	Digitoxin standard	Liquid
Quantex Digoxin	Latex-immunoturbidimetry	T3000-2275	2x9+2x6.5 mL	156	Digoxin standard	Liquid
Quantex Gentamicin	Latex-immunoturbidimetry	T3000-2281	2x12.5+2x3.8 mL	94	Gentamicin standard	Liquid
Quantex Phenobarbital	Latex-immunoturbidimetry	T3000-2278	2x14.5+2x4 mL	95	Phenobarbital standard	Liquid
Quantex Phenytoin	Latex-immunoturbidimetry	T3000-2279	2x14.5+2x4 mL	95	Phenytoin standard	Liquid
Primidone	Immunoassay	W150291	1x25+1x8 mL	106	Primidone calibrator	Liquid
Quantex Theophylline	Latex-immunoturbidimetry	T3000-2282	2x15+2x2.5 mL	95	Theophylline standard	Liquid
Quantex Valproic Acid	Latex-immunoturbidimetry	T3000-2280	2x15+2x3.8 mL	142	Valproic Acid standard	Liquid

DRUGS OF ABUSE (not available in all countries)

Test	Method	Cat. No.	Kit Configuration	Tests/Kit	Calibrators	Formulation
Amphetamines	Enzyme immunoassay	W150017	1x100+1x100 mL	952	Negative, Multi drug calibrators	Liquid
Barbiturate	Enzyme immunoassay	W150225	1x100+1x100 mL	952	Negative, Multi drug calibrators	Liquid
Benzodiazepine	Enzyme immunoassay	W150039	1x100+1x100 mL	952	Negative, Multi drug calibrators	Liquid
Cannabinoid	Enzyme immunoassay	W150185	1x100+1x100 mL	952	Cannabinoid Calibrators	Liquid
Cocaine Metabolite	Enzyme immunoassay	W150055	1x100+1x100 mL	952	Negative, Multi drug calibrators	Liquid
Cotinine	Enzyme immunoassay	W150394	1x100+1x100 mL	952	Cotinine Calibrator	Liquid
Ecstasy	Enzyme immunoassay	W15100075	1x100+1x100 mL	952	Negative, Ecstasy Calibrators	Liquid
Ethyl Alcohol	Kinetic	W150037	1x100+1x100 mL	952	Ethyl Calibrators	Liquid
Methadone	Enzyme immunoassay	W150596	1x100+1x100 mL	952	Negative, Multi drug calibrators	Liquid
Methaqualone	Enzyme immunoassay	W150514	1x100+1x100 mL	952	Negative Ecstasy Calibrators	Liquid
Opiate	Enzyme immunoassay	W150135	1x100+1x100 mL	952	Negative, Multi drug calibrators	Liquid
Phenciclidine	Enzyme immunoassay	W150160	1x100+1x100 mL	952	Negative, Multi drug calibrators	Liquid
Propoxyphene	Enzyme immunoassay	W150432	1x100+1x100 mL	952	Negative, Multi drug calibrators	Liquid