

# Avanti® J-E Centrifuge System

Compact High Performance Centrifuge that fits directly in your lab




## SPECIAL PACKAGE PRICE

### For Avanti® J-E Centrifuge System And 2 JA-Series Rotors

- Max Speed: 21,000 rpm (53,300 x g)
- Max Capacity: 4 x 1000 ml
- Maintenance Free Brushless Induction Drive
- Accel / Decel rates: 2/3
- Set Temperature Range: -10 to 40°C
- Automatic Door Interlock
- Imbalance Detector
- Dimensions : 87x76x64 cm



Small Footprint



SR Drive for Faster Accel/Decel



Extensive Rotor Library



Uses Less Energy

**This Offer is Available Only Until October 31<sup>st</sup> 2009\***

\* This offer is limited to Beckman Coulter EMEIA-Emerging Markets Area

# Avanti® J-E Centrifuge System

## Ordering Information:

PN	Description	RPM	Force
<input checked="" type="checkbox"/> 369003	Avanti JE High Performance Centrifuge	21 000rpm	53 300 x g
<b>Choose your 2 Rotors:</b>			
<input type="checkbox"/> 369687	JA-10 Fixed Angle Rotor - 6 x 500ml	10 000rpm	17 700 x g
<input type="checkbox"/> 360993	JA-12 Fixed Angle Rotor - 12 x 50ml Conical	12 000rpm	23 200 x g
<input type="checkbox"/> 339247	JA-14 Fixed Angle Rotor - 6 x 250ml	14 000rpm	30 100 x g
<input type="checkbox"/> 369691	JA-17 Fixed Angle Rotor - 14 x 50ml	17 000rpm	39 700 x g
<input type="checkbox"/> 369679	JA-18 Fixed Angle Rotor - 10 x 100ml	18 000rpm	47 900 x g
<input type="checkbox"/> 334831	JA-20 Fixed Angle Rotor - 8 x 50ml	20 000rpm	48 400 x g
<input type="checkbox"/> 342095	JA-20.1 Fixed Angle Rotor - 32 x 15ml	20 000rpm	51 500x g
<input type="checkbox"/> 334845	JA-21 Fixed Angle Rotor - 18 x 10ml	21 000rpm	50 400 x g
<input type="checkbox"/> 363055	JA-25.50 Fixed Angle Rotor - 8 x 50ml	21 000rpm	53 300 x g



**This Offer is Available Only Until October 31<sup>st</sup> 2009\***

\* This offer is limited to Beckman Coulter EMEA-Emerging Markets Area