














Iris IQ 200

Mikroskopska identifikacija delcev (sedimentov v urinu)

Kategorija	Kratika	Posnetek (primer)
Eritrociti	RBC	
Levkociti	WBC	
Skupina levkocitiv	WBCC	
Pločaste epitelne celice	SQEP	
Ostale epitelne celice	NSE	
Bakterije	BACT	
Kristali	CRYS	
Hijalini cilindri	HYA	
Neklasificirani cilindri	UNCC	
Glivice	YST	
Sperma	SPRM	
Sluz	MUC	